



INNOVATIVE BIOCIDES

Product information

Innotol CLEAN Suction System Cleaner

2500 ml REF 5011K



Characteristics

- Effectively removes biofilms and deposits in the system
- Low-foam cleaning concentrate
- Very good cleaning performance and very good material tolerance
- Prevents and effectively removes new deposits and blockages

Product type

Low-foam cleaning concentrate for the cleaning of dental suction units, incl. their drainage pipes.

Areas of application

The suction unit cleaner is a low-foam cleaning concentrate for the cleaning of dental suction units, incl. their drainage pipes. Regular use prevents the formation of deposits and blockages, and removes existing residue.

Application

The special cleaning agent is used in the event of problems as a 10% solution (100ml/1 water). The application time is 1-2 hours. If used regularly, the cleaner can be applied as a 5% solution (50ml/1 water) once or twice a week for 30 minutes. 2 litres of the ready-for-use solution is sucked up via the suction hose and the apparatus is rinsed with one litre of water after the application time. (This corresponds to 2 x ½ dosage caps in 2 litres of water).

Composition

In accordance with 648/2004/EC:

<5% amphoteric surfactants,

<5% phosphonates

Chemical-physical data

Density: 1.2 g/cm³

pH-value at 20° C: <2

Appearance: Clear,
colourless
solution

Dilution table

Working solution	5.0%	10.0%
1 litre	50 ml	100 ml
2 litres	100 ml	200 ml
3 litres	150 ml	300 ml
5 litres	250 ml	500 ml
10 litres	500 ml	1000 ml

Safety guidelines

H 319 Causes severe eye damage.

Wash thoroughly after handling. Wear protective gloves / eye protection. IN THE EVENT OF CONTACT WITH THE EYES: Rinse carefully for a few minutes with water. Where applicable remove contact lenses, if it is easy to do so. Continue to rinse. If eye irritation persists: Get medical advice/attention.



Attention

Available containers

- Plastic bottle

Contact

PRISMAN Pharma International AG
Am Stalden 16
4622 Egerkingen
Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

www.innocid.com

info@prismanpharma.com

Manufacturer

PRISMAN GmbH
Otto-Hahn-Ring 6-18
64653 Lorsch
Germany