



INNOVATIVE BIOCIDES

Product information

Innolin RAPID Spray and Wipe Disinfectant

1000 ml	REF 3010I
2500 ml	REF 3010K
5000 ml	REF 3010L
10000 ml	REF 3010N



- Unique alcohol reduced spray and wipe disinfection
- Mild and compatible with interior surfaces and users
- Comprehensive and broad spectrum of effectiveness in just 1 minute of application
- Inspected and certified according to the latest guidelines
- Licensed as class II medical device and biocide

Product type

Ready-to-use spray and wipe disinfectant based on alcohols and effective additives. Free of aldehydes and with a very pleasant and especially mild scent.

Applications

According to MPG:

Ready-for-use alcoholic fast disinfectant for the prophylactic surface disinfection of alcohol-resistant medical products for medical purposes. Not suitable for the terminal disinfection of invasive medical products.

Gem. Biocidal Products Regulation:

Ready-to-use alcoholic rapid disinfectant for the disinfection of alcohol-resistant small surfaces.

Application

Apply the rapid disinfectant to that is covers the whole surface. Then let the disinfectant take effect according to the stated application times. Only for use on alcohol-resistant materials. The disinfection solution can only be used to impregnate cloths suitable for the solution.

Composition

100g of the impregnating fluid contains:

39.0 g ethanol; 0.056 g
didecyldimethylammonium chloride

Chemical-physical data

Density: 0.938g/ml

pH-value: 6 - 9

Available containers

- Spray bottles
- Canisters

Application times

Relating to impregnating liquids

Bacteria according to EN 13727 and EN 13697 (dirty conditions):

Conc. 1 min.

Candida albicans according to EN 13624 and EN 13697 (dirty conditions):

Conc. 1 min.

TbB (tuberculosis) inactivation (*M. terra*) according to DGHM (dirty conditions):

Conc. 1 min.

Enveloped viruses (e.g., HBV, HIV, HCV) according to RKI / DVV, dirty conditions:

Conc. 1min.

Adenoviruses (EN 14476, clean conditions):

Conc. 1 min.

Noro viruses (EN 14476, clean conditions):

Conc. 1min

Rota viruses (EN 14476, clean conditions):

Conc. 1min

Relating to cloth wiping systems

Bacteria and Candida albicans according to VAH and EN16615 (dirty conditions):

Conc. 5min.

Certificates

Our products have been scientifically tested and tested according to the latest methods. We will gladly send you copies of our reports upon request.

Contact

PRISMAN Pharma International AG
Am Stalden 16
4622 Egerkingen
Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

www.innocid.com

info@prismanpharma.com

Manufacturer

PRISMAN GmbH
Otto-Hahn-Ring 6-18
64653 Lorsch
Germany

Safety Instructions

H 226 Flammable liquid and vapour.

Keep away from heat/sparks/naked flames/hot surfaces. No smoking. Keep the container tightly sealed. IN THE EVENT OF CONTACT WITH THE EYES: Rinse carefully for a long time with water. Remove any contact lenses if possible. Continue to rinse. Seal the containers and store them in a well-ventilated place. Dispose of the contents/container in accordance with local/regional/national/international regulations.

Flammable! Do not spray into naked flames! Flash point as per DIN 51755: approx. 27°C. Application max. 50ml per m². The total applied amount per room must not exceed more than 100ml per m² of room area. Keep the containers tightly sealed. Avoid eye contact.



Danger