

# Innocid<sup>®</sup>

INNOVATIVE BIOCIDES

## Product information

### Innolin WIPES AF (alcohol-free)

Disinfection Wipes

RFD 440

120 x 200 mm	REF 3021	small AF
280 x 300 mm	REF 3024	maxi AF
200 x 300 mm	REF 3028	FLOW AF



## Features

- An alcohol-free formulation for rapid disinfection even of sensitive materials
- a broad range of efficacy
- pH-neutral
- Approved as medical product and biocide (dual-use)

## Type of preparation

Alcohol-free wet wipes for cleaning and disinfecting alcohol-sensitive surfaces. It covers a broad spectrum of action against bacteria fungi and viruses.

## Areas of application

**As per German Medical Product Law MPG:** Aqueous wet wipes for prophylactic surface disinfection of non-invasive medical devices.

**As per Biocidal Regulation:** Aqueous wet wipes for the disinfection of surfaces.

## Notes on application

Open the bag with the impregnated wipes at the top and place it in the designated dispenser tub. Pull the wipe from the middle of the fleece roll and thread it through the star in the lid clip from below. Then press the lid firmly onto the tub and close tightly. The opening date, expiry date as well as the legally relevant information on the impregnating fluid must be observed and transferred to a label for external identification of the tub. After removing a wipe, the container must be closed carefully and tightly again. Carefully wipe the areas to be disinfected and allow the disinfectant to take effect. Before applying the product on large-areas of sensitive materials such as artificial leather, first test the material tolerance on a small hidden area. To prevent streaks, the treated areas should be wiped with a damp (with water) disposable cloth at regular intervals. After

## Composition

100g impregnated fluid contains:

0,24g Benzyl-C12-16-alkyldimethylammonium chloride

0,24g N,N-Didecyl-N,N-dimethylammonium chloride

0,24g C12-14-

Alkyl[(ethylphenyl)methyl]dimethylammonium chloride

## chemical-physical Data

Density: 1 g/cm<sup>3</sup>

pH value: 7,2

opening the bag, the wipes in the dispenser box can be used for 28 days. The preparation of the dispenser tub must be carried out with a suitable disinfectant as per the preparation instructions opposite. Suitable gloves (e.g. made of nitrile rubber) must be worn during application.

## Application time

Bacteria as per EN 13727 and EN 13697 (dirty conditions):

Conc. 1min.

Candida Albicans as per EN 13624 and EN 13697 (dirty conditions):

Conc. 1min.

Enveloped viruses (e.g. HBV, HIV, HCV) as per RKI/DVV, dirty conditions:

Conc. 1min.

Rotaviruses as per RKI/DVV without contamination:

Conc. 1min.

Noroviruses as per EN 14476 (clean Conditions):

Conc. 30min.

Efficacy of the impregnated wipes as per EN 16615 and VAH (dirty conditions ):

Conc. 1min.

Polyoma-SV 40 viruses (RKI/DVV with contamination):

Conc. 30min.

Impregnating fluid: Bacteria und Candida Albicans entry in VAH list:

Conc. 5 Minutes (shortest application time for surface disinfection solution in VAH list)

## Expertise

Our products have been scientifically examined and tested in accordance with the latest methods. We will be pleased to send you copies of our assessments on request.

## Contact

PRISMAN Pharma International AG

Am Stalden 16

4622 Egerkingen

Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

[www.innocid.com](http://www.innocid.com)

[info@prismanpharma.com](mailto:info@prismanpharma.com)

## Manufacturer

PRISMAN GmbH

Otto-Hahn-Ring 6-18

64653 Lorsch

Germany

## Notes on Disposal

Only put the packaging into the recycling system if it is completely emptied. Dispose of residual product in accordance with EAK 07 06 99.